


AlignShift: Bridging the Gap of Imaging Thickness in 3D Anisotropic Volumes

Jiancheng Yang*, Yi He*, Xiaoyang Huang, Jingwei Xu,
Xiaodan Ye, Guangyu Tao, Bingbing Ni.

Oct 2020



Thin-Slice (1mm)



Thick-Slice (5mm)




Axial

Sagittal

Coronal

Is there a **unified** approach to bridge the **performance gap** between thin- and thick-slice 3D **anisotropic** volumes?





TSM: Temporal Shift Module for Efficient Video Understanding. Ji Lin et al. ICCV'19.
Reinventing 2D Convolutions for 3D Images. Jiancheng Yang et al. arXiv, 1911.



Pretrained
2D Network




Pretrained
3D Network


Transfer a pretrained 2D network into a **thickness-aware** 3D network via *AlignShift*.

TSM

- Inputs are **thickness-normalized** via spatial normalization.




Original Tensor



Temporal Shift Module

AlignShift

- Inputs are not thickness-normalized.
- The converted networks are thickness-aware, i.e., the thickness is used as input.




AlignShift



3D Backbone Converted from
Truncated DenseNet-121

3D Backbone with 2D output on the **key slice**



2D Backbone	Conv2D $K \times K$	Pool2D $K \times K$	Norm2D
3D Backbone	(AlignShift+)Conv3D $1 \times K \times K$	Pool3D $1 \times K \times K$	Norm3D





Experiments: DeepLesion Benchmark



3D context modeling itself is very effective given multiple slices.

Methods	Slices	0.5	1	2	4	8	16	Avg.[0.5,1,2,4]
3DCE [26] MICCAI'18	×27	62.48	73.37	80.70	85.65	89.09	91.06	75.55
ULDor [23] ISBI'19	×1	52.86	64.8	74.84	84.38	87.17	91.8	69.22
V.Attn [24] MICCAI'19	×3	69.10	77.90	83.80	-	-	-	-
Retina. [35] MICCAI'19	×3	72.15	80.07	86.40	90.77	94.09	96.32	82.35
MVP [13] MICCAI'19	×3	70.01	78.77	84.71	89.03	-	-	80.63
MVP [13] MICCAI'19	×9	73.83	81.82	87.60	91.30	-	-	83.64
MULAN [27] MICCAI'19	×9	76.12	83.69	88.76	92.30	94.71	95.64	85.22
w/o SRL [27] MICCAI'19	×9	-	-	-	-	-	-	84.22
Ours 2.5D	×3	71.27	79.82	86.30	90.61	93.75	95.70	82.00
Ours 2.5D	×7	72.66	81.45	87.07	90.98	93.40	95.30	83.04
Ours TSM	×3	70.24	79.52	86.28	90.90	94.06	96.09	81.73
Ours TSM	×7	75.98	83.65	88.44	92.14	94.89	96.50	85.05
Ours <i>AlignShift</i>	×3	72.90	80.74	87.15	91.92	94.85	96.48	83.18
Ours <i>AlignShift</i>	×7	79.40	85.50	90.09	93.26	95.24	96.66	87.06

previous
state of the art





Experiments: DeepLesion Benchmark



AlignShift improves the performance further, which significantly outperforms previous SOTA.

Methods	Slices	0.5	1	2	4	8	16	Avg.[0.5,1,2,4]
3DCE [26] MICCAI'18	×27	62.48	73.37	80.70	85.65	89.09	91.06	75.55
ULDor [23] ISBI'19	×1	52.86	64.8	74.84	84.38	87.17	91.8	69.22
V.Attn [24] MICCAI'19	×3	69.10	77.90	83.80	-	-	-	-
Retina. [35] MICCAI'19	×3	72.15	80.07	86.40	90.77	94.09	96.32	82.35
MVP [13] MICCAI'19	×3	70.01	78.77	84.71	89.03	-	-	80.63
MVP [13] MICCAI'19	×9	73.83	81.82	87.60	91.30	-	-	83.64
MULAN [27] MICCAI'19	×9	76.12	83.69	88.76	92.30	94.71	95.64	85.22
w/o SRL [27] MICCAI'19	×9	-	-	-	-	-	-	84.22
Ours 2.5D	×3	71.27	79.82	86.30	90.61	93.75	95.70	82.00
Ours 2.5D	×7	72.66	81.45	87.07	90.98	93.40	95.30	83.04
Ours TSM	×3	70.24	79.52	86.28	90.90	94.06	96.09	81.73
Ours TSM	×7	75.98	83.65	88.44	92.14	94.89	96.50	85.05
Ours <i>AlignShift</i>	×3	72.90	80.74	87.15	91.92	94.85	96.48	83.18
Ours <i>AlignShift</i>	×7	79.40	85.50	90.09	93.26	95.24	96.66	87.06

previous
state of the art





Experiments: Analysis on Thickness

AlignShift bridges the performance gap between thin- and thick-slice volumes.

Methods	Thickness	0.5	1	2	4	8	16	Avg.[0.5,1,2,4]	diff.
2.5D ×3	All	71.27	79.82	86.30	90.61	93.75	95.70	82.00	-
	Thin	72.78	80.65	87.21	90.94	93.97	95.86	82.89	+0.89
	Thick	69.88	79.16	85.51	90.48	93.65	95.65	81.26	-0.74
2.5D ×7	All	72.66	81.45	87.07	90.98	93.40	95.30	83.04	-
	Thin	75.77	83.93	88.85	92.37	94.26	95.78	85.23	+2.19
	Thick	69.76	78.96	85.75	90.03	92.67	94.99	81.13	-1.91
TSM ×3	All	70.24	79.52	86.28	90.90	94.06	96.09	81.73	-
	Thin	73.74	81.84	87.54	92.21	94.92	96.72	83.83	+2.10
	Thick	67.03	77.37	84.98	89.95	93.16	95.44	79.83	-1.90
TSM ×7	All	75.98	83.65	88.44	92.14	94.89	96.50	85.05	-
	Thin	78.76	85.53	89.67	93.48	95.61	96.68	86.86	+1.81
	Thick	73.26	81.97	87.10	90.96	94.14	96.42	83.32	-1.73
AlignShift ×3	All	72.90	80.74	87.15	91.92	94.85	96.48	83.18	-
	Thin	73.51	81.59	87.62	92.37	94.87	96.51	83.77	+0.61
	Thick	72.85	80.18	87.10	91.94	94.91	96.54	83.02	-0.16
AlignShift ×7	All	79.40	85.50	90.09	93.26	95.24	96.66	87.06	-
	Thin	80.73	86.43	91.02	93.97	95.78	96.97	88.04	+0.98
	Thick	78.27	84.74	88.89	92.47	94.67	96.38	86.09	-0.97



Convert 2D model into 3D with a **single line of code**.

```
# m is a standard pytorch model  
m = torchvision.models.resnet18(True)
```

```
m = Converter(m,  
              alignshift_conv_cfg,  
              additional_forward_fts=['thickness'],  
              skip_first_conv=True,  
              first_conv_input_channels=1)
```

```
# after converted, m is using AlignShiftConv and capable of processing 3D volumes  
x = torch.rand(batch_size, in_channels, D, H, W)  
thickness = torch.rand(batch_size, 1)  
out = m(x, thickness)
```

Reproducible **code** on DeepLesion
with trained model **snapshots**
<https://github.com/M3DV/AlignShift/>





- We **challenge spatial normalization** as a standard pre-processing approach to solve thickness issues.
- We introduce operators to transfer any pretrained 2D network into 3D.
- We further propose a novel operator ***AlignShift*** to convert any 2D pretrained network into **thickness-aware 3D network**, which bridges the gap between thin- and thick-slice data.
- Without whistles and bells, we establish a new state of the art on DeepLesion, which surpasses prior arts by considerable margins.





Thanks for Listening



Check out my home page
for materials on this study
<https://jiancheng-yang.com/>



WeChat
jekyll4168